

# ProAir Digihaler



## ProAir Digihaler® 90 mcg albuterol

- ◆ Digital inhaler with sensors that know when the inhaler is used. It also measures how well you breathe in the medicine.
- ◆ Comes with a Mobile App. You can share information from the app with your asthma provider.



Medicine	Need a	My dose
Quick Reliever (dry powder inhaler)	No	Check Asthma Action Plan

1. Holding RespiClick upright, open red cap until you hear a click
2. Exhale away from device
3. Place mouth over mouthpiece with a tight seal (do not block vent with mouth or fingers)
4. Inhale quickly and deeply
5. Hold breath as long as you can, around 10 seconds
6. Remove RespiClick from mouth and exhale slowly away from device.
7. Place cap back on until it clicks closed
8. Steps 1-7 for each dose

**Use the camera on your phone to scan this QR code for the ProAir Digihaler Website:**



There is information about:

- ◆ How to use the ProAir Digihaler
  - ◆ Scroll down to “Download Instructions for Use”
- ◆ To save money:
  - ◆ Click “Save on ProAir Digihaler”